

Time-of-Flight Neutron Reflectometry



Suzanne te Velthuis

Materials Science Division

Instrument Scientist of POSYI reflectometer

A) Introduction:


Time of Flight Neutron Reflectometry

B) Your Experiment:

Neutron Reflectivity from Polymer Films on Si

Neutron Reflectivity

$$|\mathbf{k}| = 2\pi/\lambda$$



Θ : angle of incidence

λ : wavelength

The reflectivity of the sample is measured as a function of the scattering vector \mathbf{Q}


$$\mathbf{Q} = -\mathbf{k}_i + \mathbf{k}_f$$

$$|\mathbf{Q}| = 4\pi \sin \Theta / \lambda$$

=> two concepts for neutron reflectivity measurements:

- a) fixed wavelength + variable angle
- b) variable wavelength + fixed angle

Total Reflection at Surfaces



For neutrons (and X-rays) with wavelengths of a few Å, almost all materials have an optical index slightly smaller than 1.

=> Total reflection up to a critical angle $\Theta_{\text{crit}}(\lambda)$

Refraction index:

$$n(\lambda) = k_{z2} \text{ (inside the media)} / k_{z1} \text{ (outside)}$$

Kinetic energy of a free particle:

$$E_1 = \frac{\vec{p}^2}{2m_N}$$

Inside the media with potential V , k_{z2} is (in most cases) smaller (conservation of energy):

$$\begin{aligned} \frac{\vec{p}^2}{2m_N} + V &= E_1 \\ \Rightarrow k_{z2} &= (k_{z1}^2 - 2m_N V / \vec{p}^2)^{1/2} \end{aligned}$$

Connection to microscopic properties:

Fermi pseudo potential: $V = 2\pi \vec{p}^2 N b / m_N$

with N : number density [at/cm³]

b : coherent scattering length of the nuclei in the material [fm]


Critical angle for total reflection is reached, if $E_z = V$!

$$\Theta_{\text{crit}} = \sin^{-1} \lambda (N \cdot b / \pi)^{1/2} = \cos^{-1} n$$

or

$$Q_{\text{crit}} = 4\pi \sin \Theta / \lambda = 4(\pi N \cdot b)^{1/2}$$

Calculation of the reflectivity at a potential step




Solution of the quantum mechanic problem:

Fresnel equations

$$\text{Reflectivity} \quad R = | r |^2 = | (k_1 - k_2) / (k_1 + k_2) \exp(i2k_1 z) |^2$$

$$\text{Transmission} \quad T = | t |^2 = | 2k_1 / (k_1 + k_2) \exp(i2(k_1 - k_2)z) |^2$$

Example: Potential of a multilayer




At each interface one has to take into account:

- Refraction effects
- Multiple-scattering effects



Bragg's Law for Periodic Layered Structures

constructive interference if: $2d \sin\Theta = n \lambda$



d : double layer thickness
 Θ : angle of incidence
 n : order number (0,1,2,...)
 λ : wavelength


Reflectivity of Layered Structures



Reflectivity of Magnetic Layers

Fermi pseudo potential:



$$V = 2\pi \hbar^2 N (b_n +/- b_{mag}) / m_N$$



with b_{nuc} : nuclear scattering length [fm]
 b_{mag} : magnetic scattering length [fm]
 $(1 \mu_B/\text{Atom} \Rightarrow 2.695 \text{ fm})$
 N : number density [at/cm^3]
 m_N : neutron mass


Spin“up” neutrons see a **high** potential.
 Spin“down” neutrons see a **low** potential.

Polarized Neutron Reflectivity of Layered Magnetic Structures




Bragg's Law for Periodic Layered Structures

constructive interference if: $2d \sin\Theta = n \lambda$




d: double layer thickness
 Θ: angle of incidence
 n: order number (0,1,2,...)
 λ: wavelength



example:
 antiferromagnetic coupling
 of magnetic layers

Experiments with “Soft” Matter and Liquids



Research topics:


- morphology and thermodynamics of ultrathin polymer films
- monolayer films on liquids
- liquid/liquid and liquid/solid interfaces
- structure evolution
- diluted systems

- much more !!!

Why are (rare) Neutrons an excellent Probe for Soft Matter?




Hydrogen and Deuterium Labeling



Hydrogen

$$b_{xray} = 0.282 \times 10^{-4} \text{ \AA}$$

$$b_{neutron} = -0.374 \times 10^{-4} \text{ \AA}$$



Deuterium

$$b_{xray} = 0.282 \times 10^{-4} \text{ \AA}$$

$$b_{neutron} = 0.665 \times 10^{-4} \text{ \AA}$$


**Hydrogen and deuterium have a huge difference
in the interaction strength with neutrons !**

Examples for Hydrogen/Deuterium Contrast Variation


Neutron Scattering Length Density for Polystyrene

hydrogenated (C_8H_8): $1.4 \times 10^{-6} \text{ \AA}^{-1}$

deuterated (C_8D_8): $6.4 \times 10^{-6} \text{ \AA}^{-1}$




dPS has a much higher edge
of total reflectivity !



“Double layer” system shows two periodicities:
 a) low frequency contribution from 390 Å dPS
 b) high frequency contribution from
 690 Å total film thickness (barely visible)
 The oscillations from 390 Å dPS dominate
 due to its high potential !


The Intense Pulsed Neutron Source Facilities



Neutrons are created by a 450 MeV/14 μ A (=6.3 kW) proton beam which hits a Uranium target. The protons strike the target in 80 ns long pulses with a frequency of 30 Hz. The Uranium “boils off” neutrons in spallation and fast fission processes. The fast neutrons emerging the target are slowed down by a liquid hydrogen moderator ($T_{eff} = 32$ K).


The IPNS Neutron Reflectometer POSY II

Instrument Scientist POSYII:
Rick Goyette



Beam Line	C2
Initial Flight Path	6.2 m
Final Flight Path	1.8 m
Beam Size	50 mm x (0 - 3 mm)
Detector	Linear Position Sensitive Detector 20 cm
Choppers	none
Intensity	100 neutrons/pulse
Wave-vector Range	$0-0.25 \text{ \AA}^{-1}$
Wave-vector Resolution	$1 \times 10^{-4} \text{ \AA}^{-1}$

The Filter/Collimation System of POSY II



Soller filter:

Frame overlap filter:

=>

Neutron spectrum of POSY II: $2.5 \text{ \AA} > \lambda > 16 \text{ \AA}$

reflects $\lambda > 2.5 \text{ \AA}$ to the sample


reflects $\lambda > 16 \text{ \AA}$ out of the beam

The IPNS Neutron Reflectometer POSY I

Instrument Scientist POSYI:
Suzanne te Velthuis

Scientific Associate:
Rick Goyette

Post Doc:
Abdel Al-Smadi





Beam Line	C2
Initial Flight Path	8.3 m
Final Flight Path	0.9 m
Beam Size	(0-0.3) x 25mm
Detector	Linear Position Sensitive Detector 20 cm
Choppers	none
Intensity	40 neutrons/pulse
Wave-vector Range	$0 - 0.08 \text{ \AA}^{-1}$
Wave-vector Resolution	$2 \times 10^{-4} \text{ \AA}^{-1}$

The Filter/Collimation System of POSY I




Neutron spectrum of POSY I: $2.5 \text{ \AA} > \lambda > 14 \text{ \AA}$

Frame overlap filter:
reflects $\lambda > 14 \text{ \AA}$ out of the beam



Time-of-flight Spectrum at Detector Position



Relation between wavelength λ and time of flight t_{TOF} :

$$\lambda = \frac{h}{m_n L_{TOF}} \cdot t_{TOF}$$

with:

Planck's constant $h = 6.626 \cdot 10^{-34} \text{ Js}$

neutron mass $m_n = 1.675 \cdot 10^{-27} \text{ kg}$

distance source/detector $L_{TOF} = 7.750 \text{ m}$

$$\lambda = 2.5 \text{ \AA} \Rightarrow t_{TOF} = 4898 \text{ \mu s}$$

$$\lambda = 16 \text{ \AA} \Rightarrow t_{TOF} = 31348 \text{ \mu s}$$

new pulse starts at: 33333 μs

$$\lambda (\text{\AA}) = 5.104 \cdot 10^{-4} \cdot t_{TOF} (\mu\text{s})$$


Performing a Time-of-Flight Neutron Reflectometry Experiment



1. Step:


Measure a “Transmission Run”

(This is a measurement without a sample in the beam
in order to determine the incident wavelength spectrum.)




Looking at the Raw Data (POSYII)

Setup for a transmission run:
(top view)



Looking at the Raw Data (POSYII)




Transmission run:

plotted:
Intensity =
 $f(\text{Position on detector})$

**Straight-through
neutron beam is
in channel 236.**

**This defines the
zero point of the
reflection angle.**

INPUT ALL COMMANDS FROM TEXT WINDOW



Performing a Time-of-Flight Neutron Reflectometry Experiment



2. Step:

Insert the sample and measure a “Reflectivity Run”


(Choose the scattering angle appropriately

such that the run covers the desired Q range.)

3. Step:


**Calculate the reflectivity
of the sample**

**(Divide the “Reflectivity Run”
by the “Transmission Run”).**



Looking at the Raw Data (POSYII)


Reflectivity run:



Looking at the Raw Data

X***** []

```
20-AUG-99 INST: POSY USER: Lal           RUN/VER: 10974 /0000
17:06:34   TITLE:    2 of 2 -- dPS(390) 0.5 & 1.0 Degrees
DISPLAY: 2H-F;Y= 1, 1 :T= 1,256          PULSES: 1:103188
CURSOR:
```



INPUT ALL COMMANDS FROM TEXT WINDOW



Reflectivity run:

plotted:
Intensity =
f(Position on detector)

**Reflected neutron
beam is in channel 156.**

**The reflection angle will
be precisely calculated by
the data analysis program.**

**(The part of the neutron
beam which is not reflected
but transmitted through the
sample is absorbed by a
block of Boronnitride.)**

Performing a Time-of-Flight Neutron Reflectometry Experiment



4. Step:

Analyze the measured data to learn about your sample

(Use “Trial and Error” method but put as much previous knowledge in your models as you can !

- What are the approximate scattering length densities of your sample ?
- What are the approximate layer thicknesses ?
- etc.)



Experiment Schedule

Experiment module	Room	Assistant
Module A Sample preparation/ spin coating	Outside Chem. Lab C248	Xuesong Hu
Module B Monte Carlo simulations	A223	Suzanne te Velthuis
Module C Experiment	POSYII	R.S. Krishnan
Module D Data analysis	A223	Rick Goyette & Abdel Al-Smadi



Experiment Schedule

Time	Group 1	Group 2	Group 3	Group 4
Wed 3:30-5:30 pm	Module A	Module D	Module C	Module B
Wed 5:30-7:30 pm	Module B	Module A	Module D	Module C
Thu 2:45 -4:45 pm	Module C	Module B	Module A	Module D
Thr 4:45 - 6:45 pm	Module D	Module C	Module B	Module A

Looking at the Raw Data (POSYI)

POSYI

